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<http://openlabs.bth.se/electronics/index.php/en>

### OpenLabs Electronics Laboratory

MAIN MENU  
→ Start  
→ About  
→ Demo  
→ FAQ

OPENLABS  
BTH  
TEKNISKA HOGSKOLEN


Login  

## Welcome

Welcome to the open electronics laboratory.

Here you will find the resources needed to experiment in electronics via the internet. We have developed a system where you can make electronic experiments, right here in your browser. We supply basic equipment, such as oscilloscope, multimeter, function generator and power supply. With these and a number of electronic components you can build circuits on our virtual breadboard. None of the measurements are simulated. The circuits you build will be formed and measured on, and the real measurement results will be displayed.

Interested? Go to our [demo page](#).

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### OpenLabs Electronics Laboratory

MAIN MENU  
→ Start  
→ About  
→ Demo  
→ FAQ

OPENLABS  
BTH  
TEKNISKA HOGSKOLEN

Login  

## Demonstration

To give an idea of how the electronics laboratory is operated, a few video clips are available, demonstrating:

1. How to measure the resistance of a resistor using the multi-meter if you have problems interpreting the color code. [Show clip](#) (2.2MB)
2. How to measure U1 and U2 defined in the first circuit of an excerpt of a [lab instruction manual](#) using the oscilloscope. Please download and print the excerpt before showing this clip. [Show clip](#) (1.6MB)
3. How to perform an experiment on the operational amplifier uA741. [Show clip](#) (1.6 MB)

## Guest login

The laboratory is open for guests when its not used in regular lab sessions. Guest sessions are limited to 30 minutes.

To be able to perform the guest experiments you will have to download its [lab instruction manual](#). It describes AC experiments using passive components and experiments on an operational amplifier. When you start the experiment, you can choose from two laboratory setups at the bottom of the screen, choosing either the AC experiments or the op amp ones. This will provide all the needed components for either experiment. (No additional components need to be added).

When doing the measurements, a *virtual instructor* will check that the circuits wired will not damage components or instruments used. Please refer to the [FAQ](#) if you are not familiar with a breadboard. The oscilloscope manual can be downloaded from [manufacturers homepage](#). Please note that not all functions of the instruments are implemented.

After logging in a *virtual instrument shelf* will show up. Press Done to select the default instruments.

[Login in as a guest](#)

If you have any questions about this page or the laboratory, contact the [administrator](#).

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- Start
- About
- Demo
- FAQ

**Guest login**

In order to test out laboratory you will have to enter your email address. The information is used to track how many guests we have had.

E-Mail

If you have any questions about this page or the laboratory, contact the [administrator](#).

Clique em **Guest Course**

**OpenLabs Electronics Laboratory**

Logout 

**MAIN MENU**


- Start
- About
- Demo
- FAQ

**Welcome**

Welcome to the open electronics laboratory.

Here you will find the resources needed to experiment in electronics via the internet. We have developed a system where you can make electronic experiments, right here in your browser. We supply basic equipment, such as oscilloscope, multimeter, function generator and power supply. With these and a number of electronic components you can build circuits on our virtual breadboard. None of the measurements are simulated. The circuits you build will be formed and measured on, and the real measurement results will be displayed.

Interested? Go to our [demo page](#).

Clique em **Start Experimenting**

**OpenLabs Electronics Laboratory**

Logout 

**MAIN MENU**

- Start
- About
- Demo
- FAQ

**Guest Course**

**Start Experimenting (with flash client)**

Starting an experiment without a reservation gives you one hour of experimentation time. Experimenting without a reservation has lower priority and you will be kicked out if a reservation needs your seat.

If you have any questions about this page or the laboratory, contact the [administrator](#).

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